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REF **FB110A** 8 samples/year cod. FB, 5 mL · 4 exercises, 2 samples each.

INTENDED USE

The aim of the EQA schemes is to promote a continuous improvement of the quality of the laboratories results through retrospective comparison between labs and third party analysis of the analytical performances. Quality control samples can be used with different diagnostic kits and analyzers. For safety purposes (see Safety measures) quality control samples should be treated in the exact same manner as the patient's samples. In this way any estimated error on the control sample measurement can report a possible error on the measurement of patient samples.

PRODUCT DESCRIPTION

Control material should have chemical-physical features similar as much as possible to patients samples (commutability), the use of human matrix materials is highly recommended to ensure this requirement is met.. Samples provided are in human matrix with preservatives and stabilizers; they should be stored and treated accordingly (see: STORAGE AND STABILITY).

MATERIAL PROVIDED

Each package contains eight liquid samples, in 5 ml plastic bottles with dropper, divided according to the different Exercises of the Annual Cycle (see RESULTS FEEDBACK).

MATERIAL REQUIRED BUT NON PROVIDED

Volumetric pipettes (checked and calibrated); disposable tips; distilled water; equipment (instrument and/or kit).

STORAGE AND STABILITY

Samples stored unopened at 2-8°C in their original package are stable up to the date within which Participants should send their results back. (see: RESULTS FEEDBACK). Controls cannot be frozen.

PROCEDURE

Use samples as routine/screening samples, following the manufacturer test/instrument recommended procedure.

Dispense a sufficient control quantity in the sample container, wait at least 30 minutes before proceeding with the test which, in any case, should be run not later than 3 hours from dispensation. IMPORTANT: abide to the above mentioned schedule.

LIMITATIONS

Samples use human serum obtained from blood bags collected from Italian Blood banks and contain anti-coagulants such as: ACD (citric acid+sodium citrate+dextrose);CPD (citrate+phosphate+dextrose); CPDA (CPD+adenine). Don't use in case of evident microbiological contamination or excessive turbidity, which may cause a reduced stability. Use samples for external quality assessment only and not as standards to adjust method calibration.

SAFETY MEASURES

For in-vitro use only. The material of human origin used to manufacture samples has been tested and found non reactive with approved method for Hepatitis B surface antigen and for antibodies anti-HIV and anti-HCV. Since being negative to these tests doesn't imply that they are non infectious, to avoid personal and environmental contaminations it is recommended to treat samples as potentially infectious: use disposable gloves while handling samples and during testing; don't pipet reagent by mouth; don't smoke, eat, drink or apply cosmetics during the test. Wastes should be decontaminated and disposed of according to local regulations in force. Samples contain Sodium Azide in concentrations below 0.1% w/w. For more information on specific hazards, see related MSDS.

SENDING CONTROL SAMPLES VALUES

For each control sample, values must be sent according to the "Program calendar" available on the website. In order to send the measurements, the laboratory must log in the website www.Vanica.it, choose the EQA program and go into the "Data entry" section. The laboratory can request its credentials (username and password) by sending an e-mail at info-qc@Vanica.it.

The analytes concentrations must be reported in the units of measurement indicated in the input form; the laboratory that uses different units must convert its results using the conversion factors reported in the instructions of the diagnostic kit used or by consulting the table "Conversion factors" available on the website www.Vanica.it. We invite you to keep a print of the analyzer that carried out the measurement of the control samples.

DATA ELABORATION AND REPORTS

All the results received within the deadlines indicated in the "Calendar of programs" are included in the periodic reports. These reports are sent by e-mail according to EQA calendar to the laboratory email address, with eventual comments and supplementary analysis.

The following summary reports are prepared annually:

the summary of methods performance which contains informations about precision and accuracy of the main methods;

a summary of the laboratory's performance which contains, for each analyte, informations about laboratory accuracy and precision and the relative scores.

Laboratories can download the periodic reports in PDF format as soon as the analysis is completed by log in the website www.Vanica.it.

On the website it is also possible to download a guide to read the periodic reports, the summary of the methods performances and the summary of the laboratory performances.

In accordance with UNI CEI EN ISO 17043: 2010, the results sent by the participants are kept for 48 months.

HELP DESK

For any further information our help-desk is available from monday to friday (09:30am 5:00 pm UTC+1) at +39 010.2725628 131. Email: info-qc@Vanica.it.

CONFIDENTIAL INFORMATIONS

The communication of the EQA results between Vanica and the participants are confidential. The participants are assigned with a confidential identification code. The reports results are intended for the exclusive use of the laboratory. Any scientific publications that use these informations require the approval of Vanica.

Cod. F0262 Rev. 0 del 06/12/2022

						
SEE INSTRUCTIONS FOR USE	PRODUCT CODE	EXPIRY DATE	TEMPERATURE LIMITS	LOT NUMBER	ATTENTION	MANUFACTURER

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